

Sample Stand



Use the same sensor in the laboratory, process scale-up, and manufacturing.

This approach gives reproducible results and makes scale up from development to production easy. The analyzer is easy to move and remove from the sample stand, blender, dryer, or coater.

Specifications

Materials: All stainless steel construction

Probe fitting: 3" sanitary fitting with tri-clover clamp

Probe fitting dimension: ID 2.75" round x 2"

Base dimension: 8" square x 3 1/2" H

PART NUMBER: SS-100

Control Development, Inc., reserves the right to change specifications without notice.

Control Development, Inc.
2633 Foundation Drive, South Bend, IN 46628
Phone: (574) 288-7338 Fax: (574) 288-7339
www.controldevelopment.com
sales@controldevelopment.com

**CONTROL
DEVELOPMENT**
THE PRICE PERFORMANCE LEADER...